



tango reflex™

Setting the standard
in treatment versatility
and patient care quality

SLT FOR GLAUCOMA

POSTERIOR MEMBRANECTOMY

CAPSULOTOMY

IRIDOTOMY

Helping the world see clearly

Proprietary technologies to expand your scope of care

Every day, patients present a diverse range of eye conditions for you to evaluate and treat.

With Tango Reflex™ from Ellex, you can perform a full range of anterior and posterior YAG laser procedures.

Tango Reflex™ enables you to expand the scope of care you're able to provide, and all at the absolutely highest standard of efficacy.

A broad range of treatment options

Choose Tango Reflex™, and you can select from four treatment modalities, which comprise:

SELECTIVE LASER TRABECULOPLASTY (SLT)

POSTERIOR MEMBRANECTOMY

CAPSULOTOMY

IRIDOTOMY



tango
reflex

Offset μm



Posterior



Anterior



16

100x
0
100 P
200 P
300 P

↑ ↓ ○ ○ ○

Tango Reflex™ – treatment versatility, patient care quality

Floater & Opacities best ever visualization

Featuring Ellex's proprietary Reflex Technology™, Tango Reflex's visualization in both on-axis and off-axis modes, combined with titratable illumination, allows you to accurately visualize floaters and opacities, and to assess their position relative to the lens or retina.

SLT precision — and beyond


Tango Reflex™ incorporates Ellex's proprietary SLT technology — providing superior energy control, a sharp-edged aiming beam and the industry's fastest firing rate of three shots per second.* With an enhanced view of the trabecular meshwork, you'll be able to perform SLT procedures faster and more accurately.

YAG low energy, high efficiency

Tango Reflex's YAG mode features an Ultra Gaussian beam profile and fast rise time. That means you can perform capsulotomies and iridotomies at lower, more efficient energy levels. With less energy delivered into the eye, you'll be able to carry out capsulotomies with all types of IOLs and with significantly less risk of lens pitting*.



A proprietary dual-mode laser system, teamed with an intuitive tablet interface, allows you to instantly switch between SLT and YAG modes at the touch of a button.



The Reflex illumination mirror ensures that the laser beam is never obstructed - minimizing the risk of under-dosing or overdosing the energy and ensuring the desired therapeutic effect.

SLT; primary therapy for effective, proven glaucoma treatment

Specify Tango Reflex™, and you can deliver SLT to stimulate the eye's natural healing response to manage your glaucoma patients' IOP without the burn and scar tissue associated with argon laser trabeculoplasty (ALT).

It's a solution that, when used as a first-line option, can reduce IOP as effectively as medication, without the associated side effects or compliance problems*.

* Katz LJ, Steinmann WC, Kabir A, Molineaux J, Wizov SS, Marcellino G; SLT/Med Study Group. Selective laser trabeculoplasty versus medical therapy as initial treatment of glaucoma; a prospective, randomized trial. *J Glaucoma*. 2012;21:460-8

Stimulate safe healing

SLT delivered by Tango Reflex™ works by applying selective photothermolysis of the pigmented trabecular meshwork cells. This intervention stimulates the release of proteins and cytokines and the recruitment of macrophages to improve the circulation of aqueous humor through the meshwork and the inner wall of Schlemm's canal — all without causing thermal or coagulative damage.

Evaluate pathology

You can evaluate the result of intervention with SLT to estimate the location of the pathology. If intervention with SLT is successful, the primary obstruction region may lie within the trabecular meshwork, but if it is ineffective, then the primary obstruction may reside elsewhere.



Accurate, effective capsulotomy

A powerful solution for post-capsular fragments and opacities. Choose Tango Reflex™ and you'll secure new levels of accuracy in capsulotomy — a perfectly centered, precise capsulotomy.

It's accuracy that minimizes fringes and tags and prevents lens damage — even if the lack of a ridge makes the capsule adhere to the optic.

tango reflex.

Precision in incision

Ellex's proprietary YAG laser cavity within Tango Reflex™ delivers a four nanosecond Ultra Gaussian pulse at high peak energy and can typically achieve the industry's lowest energy optical breakdown*. This allows the energy to form a tight plasma ball, and results in less energy dispersion into surrounding tissue. This makes possible a tissue incision technique that deploys a smaller shockwave, which delivers superior treatment precision and maximizes procedure efficacy.

IOL-friendly photodisruption

Featuring an Ultra Gaussian beam profile, Tango Reflex™ focuses more energy into the center of the beam profile to deliver greater energy density — reducing the energy needed to effectively perform capsulotomy and consequently greatly reducing the risk of lens pitting*.


* Average performance only. Based on system performance testing. Data on file. Ellex Medical.

Prevent post-capsulotomy floaters

Employing TCI™ to identify capsular fragments, Tango Reflex™ can be used to vaporize broken pieces of the fragment and help prevent the common problem of sudden post-capsular floaters post treatment.



Capsulotomy, Ellex – Step 1: multifocal lens



Capsulotomy, Ellex – Step 2: multifocal lens

Specifications

	SLT Mode	YAG Mode
Laser Source	Q-switched, frequency doubled Nd:YAG	Q-switched Nd:YAG
Wavelength	green: 532 nm	infrared: 1064 nm
Energy	0.3 to 2.6 mJ per pulse, continuously variable	0.3 to 10 mJ per pulse, continuously variable
Pulse Width	3 ns	4 ns
Burst Mode	single pulse only	1, 2 and 3 pulses per burst, selectable
Spot Size	400 μm	8 μm
Cone Angle	<3 degrees	16 degrees
Offset (Anterior and Posterior)	not applicable	0, ± 100 to ± 500 μm , continuously variable
Aiming Beam	red 635 nm, adjustable intensity	green 515 nm, adjustable intensity

Common Features Specification

Repetition Rate	up to 3 Hertz
Magnification	optimized for enhanced anterior segment visualization
Cooling	Fan-cooled
Electrical Requirements	100–240 VAC, 50/60 Hz, 800 VA
Weight	31 kg, 68 lbs (laser only)
Dimensions (HxWxD)	57 x 75 x 45cm, 23 x 30 x 18 inches (laser only)
Standard Accessories	Total Solution™ tables, remote display, safety glasses, laser safety sign, dust cover
Optional Accessories	Tonometer mount, vitreolysis laser lens, SLT laser lens, capsulotomy and iridectomy laser lenses, footswitch, five-position magnification changer, beam splitter, “C” mount camera adapter, video camera adapter, co-observation tube

Specifications are subject to change without notice



tango reflex™

Find out how Tango Reflex™ will help you set new standards in SLT for glaucoma, posterior membranectomy, capsulotomy and iridotomy.

Contact us now to schedule a demonstration

Ellex Medical Pty Ltd

3 Second Avenue
Mawson Lakes, SA, 5095 AUSTRALIA
+61 8 7074 8200

Ellex Inc. (USA)

7138 Shady Oak Road
Minneapolis, MN, 55344 USA
800 824 7444

Ellex Inc. (Japan)

Harumi Center Bldg 5F, 2-5-24
Harumi Chuo-ku
Tokyo 104-0053 JAPAN

Ellex France SARL

La Chaufferie – 555 chemin
du bois
69140 Rillieux la Pape
+81 3 5859 0470

Helping the world see clearly

© Ellex Medical Pty Ltd. Tango Reflex, Reflex Technology, TCI, True Coaxial Illumination and Total Solution are trademarks of Ellex Medical Pty Ltd. Ellex is a registered trademark of Ellex Medical Pty Ltd. E&OE. International patents pending and/or granted. Tango_R 8448320EN-01

Tango Reflex™ has a CE Mark (Conformité Européenne) and US Food and Drug Administration (FDA) (510k) Market clearance for the indications of Posterior Membranectomy, Selective Laser Trabeculoplasty (SLT) Capsulotomy and Laser Iridotomy.

