



# 30XL™

| LASER SLIT LAMP |



## Optimized For Retina

The Ellex 30XL Laser Slit Lamp combines a unique 10-degree convergence angle of viewing paths with high-resolution optics to create better depth perception and an improved stereo angle for superior peripheral retina viewing.



Standard 14° stereo angle and 4mm stereo field of view

Ellex 30XL 10° stereo angle and 5.5mm stereo field of view

## Features and Options

### Total Solution™ tables



Ellex Total Solution tables are designed to meet the needs of any environment.



### Table mount

Custom mount designed to securely attach the Solitaire to a table. Adjustable tilt for easy viewing and access to controls.



### Cart

Custom designed, mobile cart provides easy access and storage for accessories.



### Transport case

Transportable case to carry and protect the Solitaire and accessories.

## Ellex 30XL Product Specifications

Eyepiece	12.5x
Magnification	5x, 8x, 12.5x, 20x, 31x
Stereoscopic Angle	10 degrees
Focal Distance, Objectives	125 mm
Light Source	12V, 30W halogen lamp
Slit width	0 - 10 mm
Aperture Diameter	0.3 mm, 1 mm, 3 mm, 5 mm, 7 mm, 10 mm
Filters	Heat Absorb (HA); Red-free (G); Cobalt Blue (B)

Specifications are subject to change without notice.

© 2007, Ellex Medical Pty Ltd. Ellex, Solitaire, Power Control, True Spot and 30XL are trademarks of Ellex Medical Pty Ltd. E&OE. PB0011B



ellex.com

### Worldwide

82 Gilbert Street  
Adelaide, SA  
5000 AUSTRALIA  
+61 8 8104 5200

### USA

7138 Shady Oak Road  
Minneapolis, MN  
55344 USA  
1 800 824 7444

### Japan

4-3-7 Miyahara 4F  
Yodogawa-ku Osaka  
532-0003 JAPAN  
+81 6 6396 2250

# Solitaire™ Solution

Premier Retina Laser for the Clinic and OR



